

NELLA™

Products For a Better Pelvic Exam

Developed with OB-GYN Providers and Patients



NELLA™
NUSPEC™

Reusable Vaginal Speculum

The Smarter Reusable Speculum



- Improved visibility with integrated sidewall retractors
- Temperature neutral—no need to warm
- Quiet operation keeps patients at ease
- Smooth—no lubricant required
- Strong, autoclavable polymer

99%

of women prefer the NuSpec
over a metal speculum¹

Full Visibility Meets Ultimate Patient Comfort


Patient: 41 y.o., BMI 30, g2p2

NuSpec



sidewall tissue retracted

Pederson



sidewall tissue encroachment

Graves



sidewall tissue encroachment

Patient and Provider Preferred.

- **Improved cervical visibility²** for patients with sidewall laxity
- **Increased comfort²** for all patients, especially adolescent, menopausal or cancer survivors
- **Minimized anxiety and discomfort²** for first-time or sensitive patients

Use NuSpec for:

General gynecology exams
Colposcopy
Biopsies
IUD insertion & removal
Other in-office procedures

Provider Needs

	NuSpec	Metal
Sidewall retention	✓	✗
Easy adjustment in use	✓	✗
Ergonomic handle	✓	✗
One-handed locking	✓	✗
No disassembly required	✓	✗

Patient Needs

	NuSpec	Metal
Smaller bill size	✓	✗
Quiet operation	✓	✗
Neutral temperature	✓	✗
Smooth material (no lubricant required)	✓	✗

¹ Ceek Women's Health (2020). NuSpec Consumer Preference Study. Internal Report: unpublished.

² Ceek Women's Health (2017). VIOLET Clinical Trial. ClinicalTrials.gov Identifier: NCT03123367. Internal Report: unpublished.

Ordering Information:

NU010 Nella NuSpec Reusable Vaginal Speculum (Box of 10)

NU001 Nella NuSpec Reusable Vaginal Speculum (Single Speculum)



NELLA™ VULIGHT®

Single-Use LED Illuminator

The Small and Mighty Speculum Light



- Shadow-free illumination
- Fits most specula, metal or plastic
- Rotates along top bill – not blocked by tools
- Reduced risk of cross contamination compared to gooseneck light

44% Brighter

The Nella VuLight has been shown to have 44% higher lux than a gooseneck exam light¹



Shadow-free:
Light where you want it.



Targeted Illumination
Nella VuLight
Single Use LED Light



Traditional Illumination
Welch Allyn 48810
Exam Light

- **Brighter illumination** compared to standard gooseneck light¹
- **Close to cervix** placement provides direct illumination even in patients with difficult anatomy, excess tissue or first-time patients
- **Securely reposition** along top bill

Use VuLight for:

General gynecology exams
LEEP procedures
Biopsies
IUD insertions/removals
Other in-office procedures
Emergency Department



Provider Needs

	VuLight	Gooseneck exam light
Focused cervical illumination ¹	✓	✗
Higher lux (brightness) ¹	✓	✗
Shadow-free ¹	✓	✗
No risk of cross-contamination	✓	✗
No need to reposition external light source	✓	✗
Portable solution	✓	✗
No provider contortion necessary	✓	✗

¹ Ceek Women's Health (2020). VuLight Comparison Verification. Internal Report: unpublished.



reddot winner 2020



Ordering Information:

VL050 Nella VuLight Single-Use
LED Illuminator (Box of 50)



NELLA™ VUSLEEVE®

Single-Use Tissue Retraction Sheath

The Only Sleeve Designed for a Speculum



- Retracts vaginal sidewalls
- Protects vaginal tissue
- Reduces need for wider specula
- Increases patient comfort¹
- Fits most specula
- Designed to grip speculum

In **80%**
of exams, VuSleeve helped
providers visualize and access the
patient's cervix better than if they
hadn't used it¹

The first-in-class, single-use tissue retraction sheath.



with VuSleeve¹
sidewall tissue retracted



without VuSleeve¹
sidewall tissue encroachment

Improve exams for all patients, including those with sidewall laxity or excess tissue

Use VuSleeve for:

- General gynecology exams
- Colposcopy
- Hysteroscopy
- Biopsies
- IUD insertion & removal
- Other in-office procedures

Provider Needs	Speculum with VuSleeve	Speculum alone
Sidewall retention ¹	✓	✗
Improved visualization & access ¹	✓	✗
Keeps procedures timely	✓	✗

Patient Needs	Speculum with VuSleeve	Speculum alone
Allows for use of smaller speculum	✓	✗
Protects sidewall	✓	✗
Reduced risk of pinching	✓	✗

¹ Ceek Women's Health (2017). VIOLET Clinical Trial. ClinicalTrials.gov Identifier: NCT03123367. Internal Report: unpublished.



reddot winner 2020



Ordering Information:

VS020 Nella VuSleeve Single-Use Retraction Sheath (Box of 20)



order@ceekwomenshealth.com or call **1-855-559-CEEK**

Product	SKU/PN	Description	UOM	Quantity	Size
Nella VuSleeve	VS020	Single-Use Universal Sheath	Box	20	OS
Nella VuLight	VL050	Single-Use LED Illuminator	Box	50	OS
Nella NuSpec	NU010	Reusable Vaginal Speculum with integrated sidewall retractors	Box	10	OS
Nella NuSpec	NU001	Reusable Vaginal Speculum with integrated sidewall retractors	Each	1	OS



Ceek Women's Health is an award-winning women led medical device company creating products for a better pelvic exam by bringing patient needs to the forefront.

ceekwomenshealth.com

1455 NW Irving Street, Suite 200 | Portland, OR 97209 USA

#ceekchange